


INSTRUCTIONAL GUIDE: REUSABLE FACE MASK WITH POCKET

Jessie Moyer

Lead Product Engineer - Trim

POCKET FACE MASK



Needed Materials


- Cotton
- 203mm of 1/4" Elastic – Qty 2)



POCKET FACE MASK



Step 1: Make paper patterns using measurements on the previous slide.



Step 2: Cut patterns out of material.

POCKET FACE MASK



Step 3: Fold and sew to notches on top. Leaving an opening in the middle.



Step 4: Press and iron seam open.

POCKET FACE MASK



Step 5: Topstitch on both sides of seam.



Step 6: Align notches on sides and iron.

POCKET FACE MASK



Step 7: Make and pin pleats.



Step 8: Iron & sew pleats into place.

POCKET FACE MASK



Step 9: Tack sew on elastic.



Step 10: Finish raw edges with binding.

POCKET FACE MASK



Mask Complete!

Disclaimer: Magna makes no warranties of any kind and does not guarantee any quantity, quality, or compliance to any standards of these masks. Magna assumes no liability for any damages or loss of any kind that might arise from the use of, misuse of, or the inability to use the masks.



DRIVING **EXCELLENCE.**
INSPIRING **INNOVATION.**